



t o i t x i ®

- INTUITIVE
- ACCURATE
- REPETITIVE

box which allows for the implantation of the entire **KAPS®** range



PL-KAP-FR-REV1 / 02 - 2017

Information

KAPS® I

L

► Indicated in cases of localized destruction of the knee joint resulting from degenerative or inflammatory involvement.

► Range of femoral and tibial components of unicompartmental knee prosthesis and associated ancillaries for mechanical replacement of the natural knee joint.

► Class III implants. France: reimbursement by the health insurance organizations according to conditions related to their inclusion on the list mentioned in Article L 165-1 of the Social Security Code.

► Before using these implants and ancillaries, healthcare professionals should carefully read the operative technique, the instructions in the package leaflet and the labeling of this product.

► 17/03/XNOV/PM/001

X·NOV

X.NO.V Medical Technology SA

Rue d'Airmont 7
2900 PORRENTUAY
SWITZERLAND

Tel. : + 41 32 544 65 60
Fax : + 41 32 544 65 79

www.xnov.com



X·NOV



KAPS®

www.xnov.com

KAPS®

1 TECHNIQUE
3 CONCEPTS
8 POSSIBILITIES



nanobiocoat®
double coating : 700µm
Titanium porous high roughness +
100µm HAP coating.

With solely one operational technique KAPS® allows the implantation of the three concepts

